

Shoprider

Rocky 6

ASK FOR THE BRAND SHOPRIDER



ROCKY 6
Thin Tyre

Available Colours



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®
www.shoprider.com.au mobility

Rocky 6

The SHOPRIDER™ **Rocky 6 Thin Tyre** is the Flagship of the SHOPRIDER range, it redefines the concept of Luxury in mobility scooters. The **Rocky 6 Thin Tyre** comes standard with 355mm (14inch) Mag wheels, Double Pillow Comfort Seat and all the features you expect from the SHOPRIDER range, fully rotating lockable seat, comfortable suspension(adjustable on rear), gas adjustable steering tiller and multi position fold down arm rests. The reliability of this impressive scooter by SHOPRIDER is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	ROCKY 6 Thin Tyre
Length	1430mm
Width	680mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Battery capacity min	2 x 60 Amp
Maximum range#	up to 50km
Wheel Size	355mm
Turning Radius	1800 mm
Weight (with battery)	146kg

FEATURES:


Double Pillow Comfort Seat
Delta Steering Tiller
Gas Adjustable Steering Tiller
Adjustable width arm rests
Adjustable rear suspension
Sliding Rotating Seat
Variable Seat Height
Battery Power indicator
Double Head Lights
Rear Lights & Indicators
Hazard Lights
Front & Rear Bump Pads
High & Low Speed
Shopping Basket
Comprehensive Warranty

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

SHOPRIDER®
www.shoprider.com.au mobility

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au

Your Authorised dealer