

BORG

Riser Recliner

“Effortless reclining
at the touch
of a button”



Upholstered in luxurious leather and matching PVC the Borg Riser Recliner Armchair is available with a single or dual motor mechanism. The dual motor will allow you to operate and move the legrest and the backrest independently. The single motor model features a spacesaver frame which only requires 6cm (2") space behind the backrest.

www.drive-medical.com.au

drive

DeVilbiss
HEALTHCARE

BORG

Riser Recliner

The Borg Range brings a contemporary style to our Restwell portfolio and these inviting chairs are truly the softest, most comfortable recliner chairs we have ever produced.



Superb Comfort, Style & Performance

PRODUCT FEATURES

- The risers have a distinct lift, rest and recline action.
- Smooth and quiet operation.
- Supreme comfort snooze position.
- Stylish design to add to any decor.
- Simple to operate handsets across the range.
- Soft and comfortable design.
- Built in storage pocket.
- Upholstered in luxurious leather with matching PVC.
- Safety key built into the hand controller.
- Battery back-up (Battery operated reset function).*
- The risers are ideal for people with limited mobility, helping them rise to their feet effortlessly.
- Safe recline distance 6cm (2") (Excludes Dual Motor).
- Features castors to aid movement.

* Batteries not included

BORG


CODE: BRC005BLKAU (Single Motor – Black)
BRC005BRAU (Single Motor – Brown)
BRC006BLKAU (Dual Motor – Black)
BRC006BRAU (Dual Motor – Brown)

Max weight:
150kg (23st)





LEATHER COLOUR OPTIONS

Black



Brown



Height:
112cm (44")



Width:
92cm (36 1/4")



Depth:
102cm (40")



Seat Depth:
55cm (21 1/2")



Seat Height:
51cm (20")



Seat Width:
51cm (20")