

# Shoprider -COUGAR

ASK FOR THE BRAND SHOPRIDER



**COUGAR POWER TILT**

**SHOPRIDER<sup>®</sup>**  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

# COUGAR

The SHOPRIDER™ **COUGAR** power tilt has all the features that come standard on the **Cougar** range of powerchair. The power tilt model features a tilting seat function which is easily operated via the joy stick controller. giving the operator a maximum recline of 45 degrees for greater user comfort. This impressive powerchair will provide smooth operation with its mid wheel drive and powerful motors. It's ideal for both inside and outside use on even ground and can travel long distances between charges.

## SPECIFICATIONS

## G424MT

Length	1090mm
Width	620mm
Height	1240mm
Battery Capacity (Minimum Amp)	2 x 50Amp
#Maximum Range	up to 40km
Maximum Load Single User*	136Kgs
Charger on board	5Amp
Motor Size	0.80 Hp
*Max Climbing Angle	9 Degrees
Turning Radius	600mm
Wheel Size:	
- Front Caster	2 x 125mm
- Middle (drive)	2 x 250mm
- Rear caster	2 x 100mm
Maximum Seat Recline	45 degree

## FEATURES:

Large Captain's Seat  
Electric power tilt seat  
Adjustable arm rests  
Swing away leg rest  
Sleek body design  
Color change panels  
Smooth operation  
Free wheel device  
Variable seat height  
Controller interchangeable  
both arms  
Off Board charger


Driver weight may exceed weight of unit. Reduce speed when turning

## IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au