

# Shoprider Crossover

ASK FOR THE BRAND SHOPRIDER



**GK10**

Available Colours



**SHOPRIDER**<sup>®</sup>  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## Crossover

**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

The SHOPRIDER™ **Crossover** is a fully portable model that offers exceptional functionality for each user. The Crossover has been designed to offer each user the flexibility of a scooter that can be easily pulled apart for transport or used just as a normal mid size scooter. The Crossover is loaded with features like 265mm wheels, mid back fully rotating captain seat, adjustable steering tiller, multi position lift up armrests & 36amp batteries.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

Length	1190mm
Width	585mm
Height	960mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.2Hp
Maximum Speed	9kmph
Battery Capacity (Minimum Amp)	2 x 36Amp
Maximum Range #	up to 27km
Wheel Size	265mm
Turning Radius	1330mm
Weight (with battery)	75kg

### GK10

### FEATURES

Mid Back Captains Seat  
 Loop Steering Tiller  
 Adjustable Steering Tiller  
 Adjustable width arm rests  
 Sliding Rotating Seat  
 Variable Seat Height  
 Battery Power indicator  
 Head Lights  
 Front & Rear Bump Pads  
 Pull Apart Design for Transport  
 Shopping Basket  
**Comprehensive Warranty**


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
 Or visit our Web Site at www.shoprider.com.au