

KING COBRA Scooter

“Powerful
performance and
dynamic design”



Following the success of our Cobra scooter we've developed the new King Cobra. This offers an evolved dynamic design, generous spaciousness, a 50km range and a maximum weight limit of 200kg.

www.drive-medical.com.au

drive

DeVilbiss
HEALTHCARE

KING COBRA Scooter

The King Cobra is a revolutionary scooter in terms of design, aesthetics and styling. It has been aerodynamically designed to give an outstandingly stable drive and maximise battery consumption. It's a remarkable scooter that makes a huge statement and sets a new benchmark with its striking dynamic design.

Superb Comfort, Style & Performance...

- Top speed of 13kmph (8mph)*
- Maximum range of 50km (32 miles) on full battery charge*
- Stylish & distinctive design
- Comfortable rotating seat with sliding and reclining facility
- Height adjustable headrest and angle adjustable armrests to suit user requirements
- Electromagnetic brake system and emergency cable brakes.
- All round suspension ensures a smooth ride and optimal handling
- Generous tiller adjustment to suit user requirements and assist easy transportation
- User friendly controls with a modern digital LCD display
- Automatic speed reduction when steering for added safety
- Twin USB power socket
- Fully functioning indicators and twin adjustable rear view mirrors for extra visibility
- High level indicators for added safety
- Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark
- Stylish 'Blade' alloy wheels with black pneumatic tyres for a smoother ride
- Delta bars as standard
- Bottle/cup holder included
- Rear anti tip wheels as standard
- Anti roll back system safety device prevents scooter from rolling backwards on hills
- Freewheel facility allows movement of scooter without switching the motor on
- Optional rear pannier box
- Recommended MK Gel batteries: 2 x 12V 96Ah (Not included)
- Weight capacity 200kg (32st)

* Subject to conditions of use

CODE: KCOBRARDAU

DESCRIPTION: King Cobra in Red



Length:
170cm (67")



Height:
140cm (54")



Width:
73cm (29")



Max. Gradient:
10°



Weight:
149kg (328lbs)



Max. Weight
Capacity:
200kg (32st)