




- LIVE YOUR BEST'-







Drive Wheels
Front Casters
Rear Anti-tips
Maximum Speed
Ground Clearance
Turning Radius
Length

Width Seat Type Batteries Range

Charger Electronics Weight Capacity

Base Weight

Standard Seat Weight

Battery Pack Weight Suspension

Heaviest Piece (rear)

Warranty

20.32 cm (8") (solid)

12.7 cm (5") (solid)

7.6 cm (3") (solid)

5.65 km/h (3.5 mph)

3.49 cm (1.37")

68.26 cm (26.87") with foot platform

82.55 (32.5") cm with foot platform

48.26 cm (19")

Foldable, compact, lightweight

(2) 12V, 12 amp

Up to 12.87 km (8 miles)

Off-board, 2 amp

40 amp, PG GC 2 (non-programmable)

113 kg (250 lbs.)

25.86 kg (57 lbs.) (front & rear)

11.34 (25 lbs.) kg

9 kg (20 lbs.)

Limited

15.65 kg (34.5 lbs.)

2-year limited

## **Performance**

**Portability** 

 Top speed of 5.65 km/h (3.5 mph) and enough power to get you where you want to go

• Easy-to-remove battery pack

· Fits in most vehicle trunks

manageable pieces

· Conveniently disassembles into four

 Twin motors for traction and precision handling

## Colours





Red

Blue

## Manoeuvrability

· Compact design for a tight turning radius

20-24 Apollo Drive, Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 6611

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

This product is not coded by DME PDAC or does not meet criteria: K0899.



