

# SERENA

## Standard Back

“Washable  
head and armrest  
covers”



The spacious Serena Riser Recliner is ideal for most users, providing a high-level-of-comfort. It features a single-motor and early-rise ottoman to enable the user to raise their legs without fully reclining the chair.

[www.drive-medical.com.au](http://www.drive-medical.com.au)

**drive**

**DeVilbiss**  
HEALTHCARE

# SERENA

## Standard Back

The Serena Riser features a contemporary design with padded upholstery and washable head and armrest cushions for added protection. The riser has a battery back-up system which allows the chair to be raised and lowered in the event of a power failure.

### Superb Comfort, Style & Performance...

- Early rise legrest before the chair is fully reclined.
- One motor activates a distinct lift, rest and recline action.
- It will rise and tilt users to their feet.
- Excellent snooze position for supreme comfort and relaxation.
- Quiet motors provide a smooth and discreet operation.
- Ergonomically contoured seat and back and padded upholstery for added comfort and support.
- Fully padded chaise for increased comfort.
- Two convenient side pockets for easy storage.
- Protective, washable head and armrest covers.
- Quality assured, meeting all relevant fire retardancy standards.
- Built-in fail-safe battery back-up system, which enables the chair to be raised and lowered in the event of a power failure\*.
- Comes with a simple-to-operate handset.
- Stitched handset holder to prevent loss or damage to the handset.
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility (e.g. those with MS, Arthritis, Rheumatism and circulation problems).
- Available in 2 colour options: Forest and Mudbrick (Micro Suede fabric).
- 2-year warranty on wooden frame and electrical components.
- Maximum weight capacity 130kg (20st/285lb).

\* Batteries included.


### SERENA STANDARD BACK

**CODE:** 920SFAU (Forest Micro Suede)  
920SMAU (Mudbrick Micro Suede)

**Max Weight Capacity:** 130kg (20st/285lb)



Height:  
111cm (43.5")



Seat Height:  
47cm (18.5")



Width:  
88cm (34.5")



Seat Width:  
54cm (21.5")




Depth:  
94cm (37")




Seat Depth:  
53cm (21")

#### COLOUR OPTIONS



Forest



Mudbrick