# VENOM Scooter Range

## "Compact in size big on features"



The Venom has many new features which are only normally found on larger scooters. It's the first Drive mini scooter to include Delta Bars as standard, it also features the next generation Drive splitting mechanism.

The Venom is ideal for use indoors and outdoors where the terrain is forgiving (eg. smooth paths, shopping centres, etc). The Venom has been developed to ensure a focus on Comfort, Style, Performance and Reliability.





www.drive-medical.com.au





## VENOM Scooter Range



#### FOUR WHEEL MINI SCOOTER (20AH)

#### COLOUR OPTIONS: MS009BLAU-1 - Blue MS009RDAU-1 - Red

Model MS009	<b>Length</b> 108cm (42.5")
Width	Height
48cm (19")	91cm (36")
Max Gradient	Weight
8°	48kg (106lb)
Max Weight	Max Range
Capacity	20km
130kg (20st)	(13 miles)

- Top speed of 6kph (4mph)\*
- Maximum range of 20km on full battery charge\*
- Comfortable lightweight seat which swivels 360°
- · On board and in-house charging as standard
- Automatic electromagnetic brake system
- Tool free height adjustable seat
- Delta Bars as standard
- Easy to change body panels
- Easy tiller adjustment
- Black puncture proof non-marking tyres
- Width adjustable armrests
- Front removable storage basket as standard
- Rear anti tip wheels as standard

### FOUR WHEEL MINI SCOOTER (12AH)

#### 1 OOK WITEEL WIIN SCOOTER (IZAIT)





Model	Length
MS008	108cm (42.5")
<b>Width</b>	Height
48cm (19")	91cm (36")
Max Gradient	Weight
8°	44kg (97lb)
Max Weight	Max Range
Capacity	16km
130kg (20st)	(10 miles)

- Top speed of 6kph (4mph)\*
- Maximum range of 16km on full battery charge\*
- Comfortable lightweight seat which swivels 360°
- Removable battery pack with easy to carry handle and locking mechanism
- On board and in-house charging as standard
- Automatic electromagnetic brake system
- Easy tiller adjustment
- Removable seat, fold down tiller and splits easily for storage and transportation
- Stylish silver wheel hubs
- Black puncture proof non-marking tyres
- Rear anti tip wheels as standard

### THREE WHEEL MINI SCOOTER



#### COLOUR OPTIONS: MS003PBAU-01 - Blue MS003RDAU-01 - Red

Model	Length
MS003	108cm (42.5")
Width	Height
48cm (19")	91cm (36")
Max Gradient	Weight
8°	44kg (97lb)
Max Weight	Max Range
Capacity	16km
130kg (20st)	(10 miles)

- Maximum range of 16km on full battery charge\*
- Top speed of 6kph (4mph)\*
- Comfortable lightweight seat which swivels 360°
- Removable battery pack with easy to carry handle and locking mechanism
- On board and in-house charging as standard
- Automatic electromagnetic brake system
- Removable seat, fold down tiller and splits easily for storage and transportation
- Stylish silver wheel hubs
- Black puncture proof non-marking tyres
- Width adjustable armrests
- $\ensuremath{^{*}}$  Depends on user mass, terrain and battery condition

www.drive-medical.com.au

Leaflet code: DDVH\_036\_Venom\_Range



<sup>\*</sup> Depends on user mass, terrain and battery condition.

 $<sup>\</sup>ensuremath{^*}$  Depends on user mass, terrain and battery condition.