

Shoprider

889AE **ROCKY 8**

ASK FOR THE BRAND SHOPRIDER



NEW MODEL!
MARCH 2020

Available Colours

**BURGUNDY
RED**

**COBALT
BLUE
(NEW COLOUR!)**

SHOPRIDER®
shoprider.com.au *mobility*

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®
shoprider.com.au *mobility*

889AE **ROCKY 8**

The SHOPRIDER™ **ROCKY 8** is the latest scooter from SHOPRIDER™ and it is loaded with great features. The ROCKY 8 has new and updated features including large low profile tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller, fold down arm rests, powerful LED lights, USB outlet, new body covers and modern electronic dash. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

ROCKY 8

Length	1590mm
Width	730mm
Height	1285mm
*Safe Climbing Angle	12 degrees
*Max Climbing Angle	12 degrees
Maximum Load Single User*	225kgs
Charger	8Amp
Motor Size - 4Pole	2Hp
Maximum Speed	10kmph
Battery Capacity (Minimum Amp)	75Amp x 2 (Included)
Maximum Range #	Up to 54km
Wheel Size	330mm
Turning Radius	1860mm
Weight (with battery)	147kg

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.

WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

NEW / UPDATED FEATURES

- Large Captain's Seat
- Adjustable Steering Tiller
- Adjustable Width Arm Rests
- Sliding / Rotating Seat
- Front Shopping Basket
- Horn
- USB Power Outlet
- LED Front Headlights
- Indicators / Hazard Lights
- High / Low Speed Switch
- Anti Tip Wheels
- Front + Rear Suspension
- Comprehensive Warranty**


LCD Display shows:

- Speed
- Battery Level
- Right / Left Turn Signal
- Hazard Warning Lights
- High / Low Speed
- Temperature
- Odometer



ROCKY 8

*(Single user)



SHOPRIDER®
shoprider.com.au *mobility*

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au

Your Authorised dealer