

Shoprider

GK9 / 3

ASK FOR THE BRAND **SHOPRIDER**



GK9/3 PORTABLE SCOOTER

SHOPRIDER®
www.shoprider.com.au
mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®

www.shoprider.com.au

mobility

GK9 / 3

The SHOPRIDER™ **GK9 / 3** is designed to be fully portable and offers the user exceptional functionality. This scooter offer many of the features found on larger scooters with the added benefit of its easy pull apart design. Features include fully rotating lockable seat, adjustable tiller adjustable height armrests, fold up armrests and 21amp battery capacity. The GK7 / 3 is highly manouvable, combine this with the legendary Shoprider reliability it will offer years of quality service and reliability.

| SPECIFICATIONS | GK9/3 | FEATURES |
|---------------------------|------------|-----------------------------|
| Length | 990mm | Mid Back Seat |
| Width | 560mm | Adjustable Steering Tiller |
| Height | 880mm | Adjustable Width Arm Rests |
| *Safe Climbing Angle | 9 Degrees | Fold up Arm Rests |
| *Max Climbing Angle | 9 Degrees | Rotating Seat |
| *Maximum Load Single user | 136Kgs | Front Shopping Basket |
| Charger (Del with min) | 3Amp | Variable seat height |
| Motor Size | 0.75Hp | Horn |
| Maximum Speed | 6Kph | Battery Level Warning Light |
| Battery capacity min | 2 x 21Amp | Battery Level Gauge |
| Maximum range # | up to 20km | Split Battery Box |
| Wheel Size | 220mm | Solid Tyres |
| Heaviest Part | 16.5kg | Anti Tips Wheels |
| Weight (w/o Batteries) | 50.0kg | Multi colour change panels |

Comprehensive Warranty


Shoprider™ is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

WARNING: Do not exceed 10kph on public roads or Paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER®

www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au