



Interceptor S2



www.meritsaustralia.com



Interceptor S2



- Brand new sporty streamline design
- Fully adjustable suspension system that delivers a smooth and comfortable ride
- Modern LED lighting system delivers a brilliant light with less power consumption
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Built-in half speed switch which operates automatically when turning to provide extra safety
- With Merits 4PT (4 pole technology) and large controller delivering much more power for climbing steep driveways and hills
- With new easy to access 12 volt socket







Modern LED lighting system delivers brilliant light with less power consumption

Dimensions cm/in









More space for storage





Digital LCD dashboard display Handy 12 volt socket

Technical Specs

Max Speed	14/8.9mph	
Range	40-60km/26-38mi	
Turn Radius	180cm/70.9"	
Ground clearance	13cm/5.1"	
Motor	DC 24V, 650W (Rating)/ 2200W (Max)	
Controller	PG S-Drive 200A	
Battery	12V/60-100A/H x 2pcs	
Charger	8A off-board	
Gradient	10°	
Front Wheel	14" (3.50-8) pneumatic	
Rear Wheel	16" (4.00-10) pneumatic	
Brake	Intelligent, regenerative, electromagnetic brakes	

Australian Standards Compliant AS/NZS 3695.2:2013

Overall Dimension: A x B x C 154 x 72 x 145.5cm/61 x 29 x 57.3"

Seat Dimension: D x E x F 51 x 46 x 63cm/20 x 18 x 25"

Seat Height from Ground: G 67.1cm/26.4"

Seat Height from Deck: H 48.6cm/19"

Wheel Base: I 113cm/45"

Weight kg/lb

Base	Battery/ea	Seat	Total with Battery	Capacity
84/185	25/55	18/40	165/364	204/450



www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.