

Shoprider

QT8-3

ASK FOR THE BRAND **SHOPRIDER**



Available Colours



SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER[®]
www.shoprider.com.au *mobility*

QT8-3

The SHOPRIDER™ **QT8-3** is a versatile three wheel portable scooter that offers the performance and features that you expect from a SHOPRIDER™ scooter. This portable scooter is easily disassembled for transport with NO requirement for tools. The **QT8-3** has been designed to offer a number of features like LED battery power display, larger leg room, super maneuverable steering and improved comfortable swivel seat. The reliability of this impressive portable scooter by SHOPRIDER is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

QT8-3

Length	960mm (w-basket 080mm)
Width	560mm
Height	985mm
Maximum Load Single User*	113Kgs
*Max Climbing Angle	6 Degrees
Motor Size	0.3Hp
Charger	2 Amp
Maximum Range #	up to 18km
Maximum Speed	6Kph
Battery Capacity (Minimum Amp)	12 Amp
Wheel Size	200mm
Turning Radius	830mm
Heaviest part (front chassis)	12kg
Weight (with battery)	41kg

FEATURES

Mid Back Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Rotating Seat
Folding Arm Rests
Variable Seat height
Front Shopping Basket
Horn
Flat Free Tyres
LED Battery Level indicator
Great Turning Circle
Easy to Dismantle & Reassemble
Anti Tips Wheels
Seperate Charger
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

QT8-3



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au