

Shoprider Crossover

ASK FOR THE BRAND **SHOPRIDER**



GK10

Available Colours



SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

Crossover

SHOPRIDER[®]
www.shoprider.com.au
mobility

The SHOPRIDER™ *Crossover* is a fully portable model that offers exceptional functionality for each user. The Crossover has been designed to offer each user the flexibility of a scooter that can be easily pulled apart for transport or used just as a normal mid size scooter. The Crossover is loaded with features like 265mm wheels, mid back fully rotating captain seat, adjustable steering tiller, multi position lift up armrests & 36amp batteries.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

	GK10
Length	1190mm
Width	585mm
Height	960mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.2Hp
Maximum Speed	9kmph
Battery Capacity (Minimum Amp)	2 x 36Amp
Maximum Range #	up to 27km
Wheel Size	265mm
Turning Radius	1330mm
Weight (With Battery)	78kg
Weight (Battery Each)	11kg (22kg Total)
Weight (Chair)	16kg
Weight (Rear Transaxle)	17kg
Weight (Main Chassis)	23kg

FEATURES

- Mid Back Captains Seat
- Loop Steering Tiller
- Adjustable Steering Tiller
- Adjustable width arm rests
- Sliding Rotating Seat
- Variable Seat Height
- Battery Power indicator
- Head Lights
- Front & Rear Bump Pads
- Pull Apart Design for Transport
- Shopping Basket
- Comprehensive Warranty**


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Your Authorised dealer

For your nearest Distributor call 1800 672 633
 Or visit our Web Site at www.shoprider.com.au